

ILEOSTOMY PRODUCTS



All the products featured in this brochure are available through Salts Medilink® or your usual supplier of ostomy products.

For more information and advice about products featured in this brochure, Stoma Solutions or the Salts Medilink® 24* Hour Home Delivery Service, please contact:

FREEPHONE 0800 626388 (UK) or **FREEPHONE 1-800 408508** (Ireland)
or log on to our website at www.salts.co.uk/medilink



100% designed and manufactured in the UK

For free samples please call

FREEPHONE 0800 626388 (UK) or

FREEPHONE 1-800 408508 (Ireland)

or log on to our website at

www.saltsstomacare.co.uk

SALTS
HEALTHCARE

© Registered trademarks of Salts Healthcare Ltd. © Salts Healthcare Ltd 2010. Products and trade marks of Salts Healthcare Limited are protected by UK and foreign patents, registered designs and trade marks. For further details, please visit www.salts.co.uk *If required, except Highlands and Islands and subject to availability. RM763886 04/10/3K



A complete guide to our comprehensive range of ileostomy products for the one- and two-piece user



Validating Salts research
into healthy stoma skin

SALTS
HEALTHCARE

Excellence in stoma care

A range of clipless, drainable pouches for the one-piece user, incorporating a new anatomical shape; thicker, silky soft covering; oval-shaped wafer; and quieter films – our most comfortable drainable pouch ever.



DRAINABLE

SECURE

New oval-shaped wafer with thumb tab, for easy removal. All starter holes cut up to 65mm. (Large transparent pouch cuts up to 70mm.)

LARGE OVAL WAFER

DISCREET
Reduced headroom – top of pouch does not fold over.

NEW SHAPE
Slightly anatomical shape, for easy draining and increased comfort.

SOFT

New thicker, silky soft covering with a new softer weld along edge of pouch and a new gently rounded profile to outlet.

QUIET
New quieter inner layers with reduced rustle, for increased confidence.

WIDE CHOICE

Available in a wide range of sizes – either starter hole or pre-cut – with the choice of opaque, transparent or transparent with overlapping fabric.



LARGE POUCHES						
SIZE	STARTER HOLE	PRE-CUT	ORDER CODES			PACK SIZE
			Transparent	Trans/Overlap	Opaque	
13mm*	•	–	–	DLTO13	–	30
13–65mm	•	–	–	–	CFDL13	30
13–70mm	•	–	CFDLT13	–	–	30

STANDARD POUCHES						
SIZE	STARTER HOLE	PRE-CUT	ORDER CODES			PACK SIZE
			Transparent	Trans/Overlap	Opaque	
13–65mm	•	–	–	CFDT13	CFD13	30
25mm	–	•	–	–	CFD25	30
28mm	–	•	–	–	CFD28	30
32mm	–	•	–	–	CFD32	30
35mm	–	•	–	–	CFD35	30
38mm	–	•	–	–	CFD38	30
13mm*	•	–	–	DTO13	–	30

SMALL POUCHES						
SIZE	STARTER HOLE	PRE-CUT	ORDER CODES			PACK SIZE
			Transparent	Trans/Overlap	Opaque	
13–65mm	•	–	–	–	CFDS13	30

* Oval wafer pouches can be cut up to 110mm wide and 70mm deep.

A drainable convex pouch with even more flexibility and softness. Confidence Convex Supersoft incorporates a soft, flexible wafer with a gentle profile – for increased comfort.



DRAINABLE

GREATER FLEXIBILITY

The wafer has an incredibly soft and flexible foam backing, producing a better fit over the contours of the skin, while retaining full and unrestricted movement.



ADDED COMFORT

Softer, more flexible convex insert, for added comfort.

MAXIMUM SECURITY

Soft, integrated 'hook and loop' opening/closing system, for maximum security.

EASIER EMPTYING AND CLEANING

New soft and flexible offset strips allow easier draining and cleaning of the outlet.



EXTRA SECURITY

New hydrocolloid wafer with stickier adhesive and a shallower profile, for extra comfort.

BELT LOOPS

Belt loops, for extra security.

SOFT

Soft beige cover, for extra comfort.

NEW SHAPE

New anatomical shape pouch.

DISCREET

Closure can be hidden under soft cover, for extra comfort and discretion.

LARGE DRAINABLE POUCHES

SIZE	STARTER HOLE	PRE-CUT	ORDER CODES Trans/Overlap	PACK SIZE
13–25mm	•	–	CDSSL1325	10
13–38mm	•	–	CDSSL1338	10
13–52mm	•	–	CDSSL1352	10

STANDARD DRAINABLE POUCHES

SIZE	STARTER HOLE	PRE-CUT	ORDER CODES Trans/Overlap	PACK SIZE
13–25mm	•	–	CDSS1325	10
13–38mm	•	–	CDSS1338	10
13–52mm	•	–	CDSS1352	10
21mm	–	•	CDSS21	10
25mm	–	•	CDSS25	10
28mm	–	•	CDSS28	10
32mm	–	•	CDSS32	10
35mm	–	•	CDSS35	10
38mm	–	•	CDSS38	10
41mm	–	•	CDSS41	10

SMALL DRAINABLE POUCHES

SIZE	STARTER HOLE	PRE-CUT	ORDER CODES Trans/Overlap	PACK SIZE
13–25mm	•	–	CDSSS1325	10
13–38mm	•	–	CDSSS1338	10
13–52mm	•	–	CDSSS1352	10

Note: Confidence Convex Supersoft should only be used after prior assessment by a stoma care specialist.

Our very popular range of cliplless pouches for the one-piece user, incorporating a soft and comfortable beige cover designed for easy cleaning, maximum comfort and security.

UNIQUE FILTRATION SYSTEM

High-performance filter and help alleviate 'ballooning'. T repelling cover to prevent od moisture out.



SECURITY

Large diamond-shaped wafer eases adhesion and wear e. Suitable for even the most sensitive skin.

NO CLIP REQUIRED

Soft and integral, easy-to-use 'hook and loop' opening/ closing system, for maximum comfort and security.



Closure hidden under soft cover, for extra comfort and discretion.

EASY CLEANING

Shorter inner layers and unique flexible strips that form a spout, for easy draining and cleaning.

STANDARD POUCHES

SIZE	STARTER HOLE	PRE-CUT	ORDER CODES			PACK SIZE
			Transparent	Trans/Overlap	Opaque	
13mm-52mm	•	-	-	SST13	-	30
25mm	-	•	-	SST25	-	30
32mm	-	•	-	SST32	-	30

A range of drainable pouches for the one-piece user, incorporating a soft and comfortable beige cover with an advanced filter designed to eliminate odour and moisture, and help reduce pouch 'ballooning'. Pouches are available in a choice of opaque fabric, or transparent with overlapping fabric.

SECURITY

A generous diamond-shaped wafer, to increase adhesion and wear time. Suitable for even the most sensitive skin. Its shape and backing make it exceptionally easy to use, comfortable and very secure.

EASY CLEANING

Shorter inner layers to the pouch opening, for easier cleaning.



UNIQUE FILTRATION SYSTEM

An innovative combination of high-performance filter and multiple layers of film allows air to flow at more than double the rate of most available drainable pouches, to alleviate the problems associated with 'ballooning'.

The filter has a water-repelling cover to prevent odour escaping and keep moisture out.

STANDARD POUCHES

SIZE	STARTER HOLE	PRE-CUT	ORDER CODES			PACK SIZE
			Transparent	Trans/Overlap	Opaque	
13-52mm	•	-	-	GDTF13	GDF13	30 + 2 clips

SMALL POUCHES

SIZE	STARTER HOLE	PRE-CUT	ORDER CODES			PACK SIZE
			Transparent	Trans/Overlap	Opaque	
13-52mm	•	-	-	-	GDSF13	30 + 2 clips

Harmony[®] Duo

All the advantages of a one-piece system, such as an ultra slim profile and flexible wafer, with all the convenience of a two-piece system: its unique 'docking system' and skin-friendly flange makes pouch changes so much easier.



DRAINABLE

SKIN FRIENDLY

The soft, flexible wafer can stay in place for longer than a one-piece, reducing skin trauma caused by frequent changes.

SOFT AND FLEXIBLE

The low-profile flanges are soft and flexible, designed to move as you do, and are available in a choice of sizes, and either standard or with hydrocolloid surround.

EASY APPLICATION

A locator on the wafer matches up with a notch on the pouch for simple, accurate positioning every time.

LIGHTWEIGHT COVER

Feels soft on the skin, and is discreet under the majority of clothing.



ANATOMICALLY SHAPED

The pouches are anatomically shaped to fit with the contours of your body, for extra comfort and freedom of movement.

EASIER EMPTYING AND CLEANING

New soft and flexible offset strips allow easier draining and cleaning of the outlet.

MAXIMUM SECURITY

Soft, integrated 'hook and loop' opening and closing system, for maximum security.

DISCREET

Closure can be securely held in place under soft cover by hidden hook material, for extra comfort and discretion.



STANDARD FLANGES

SIZE	ORDER CODES	PACK SIZE
13-32mm	HD1332	10
13-50mm	HD1350	10
13-70mm	HD1370	10

FLEXIBLE FLANGES (with hydrocolloid surround)

SIZE	ORDER CODES	PACK SIZE
13-32mm	HDF1332	10
13-50mm	HDF1350	10
13-70mm	HDF1370	10

LARGE DRAINABLE POUCHES

SIZE	ORDER CODES		PACK SIZE
	Transparent	Trans/Overlap	
13-32mm	-	HDDL1332	30
13-50mm	HDDL1350	HDDL1350	30
13-70mm	HDDL1370	HDDL1370	30

STANDARD DRAINABLE POUCHES

SIZE	ORDER CODES Trans/Overlap	PACK SIZE
13-32mm	HDD1332	30
13-50mm	HDD1350	30
13-70mm	HDD1370	30

SMALL DRAINABLE POUCHES

SIZE	ORDER CODES Trans/Overlap	PACK SIZE
13-32mm	HDDS1332	30
13-50mm	HDDS1350	30
13-70mm	HDDS1370	30

Confidence Comfort Paediatric has all the benefits of our adult product and is really kind to young, delicate skin. Confidence Comfort Paediatric is even easier to open, drain and clean. The flange is oval-shaped, for added comfort, and can be positioned so that it drains to the side – perfect for a good night's sleep.

OVAL-SHAPED FLANGE

allows drainage to the side if the child is lying down.

SOFT AND GENTLE

with a silky soft covering, softer edges and outlet profile.

NO HOLE WAFER OPTION

for flexible positioning. Hole does not have to be cut in centre of wafer.

EASY TO USE – soft, flexible outlet is comfortable and very easy to empty and clean.



KINDER TO DELICATE YOUNG SKIN

– proven, skin-friendly, non-allergenic hydrocolloid.

CHOICE OF COVER

– transparent or transparent with overlap.



PAEDIATRIC POUCHES			
SIZE	ORDER CODES		PACK SIZE
	Transparent	Trans/Overlap	
8–60mm Starter	CPT8	CP8	30
NO HOLE	CPTNH	CPNH	30

FRESHAIRE

FreshAire is effective and easy to use. During pouch changes, FreshAire is simply sprayed once or twice into the air around the pouch.



PERI-PREP[®] SENSITIVE

Leakage and the trauma caused by the constant changing of a pouch can lead to sore skin. Peri-Prep Sensitive is specifically designed to leave a gentle protective film on the skin immediately around the stoma.



SECUPLAST[®] HYDRO

Skin-friendly hydrocolloid strips to provide extra security around the edge of the flange.



WIPEAWAY

WipeAway sachets are used to remove pouches and adhesive residue and to help keep the skin around the stoma healthy. The natural ingredients help to soothe skin trauma associated with adhesive removal, leaving a hint of fresh lavender.



FRESHAIRE ODOUR ELIMINATING SPRAY			SECUPLAST [®] HYDRO		
SIZE	ORDER CODE	PACK SIZE	SIZE	ORDER CODE	PACK SIZE
50ml	FA1	1	One size fits all	SPH1	30
PERI-PREP [®] SENSITIVE			WIPEAWAY ADHESIVE REMOVER		
SIZE	ORDER CODE	PACK SIZE	SIZE	ORDER CODE	PACK SIZE
–	PPS1	30 sachets	–	WA1	30 sachets

DERMACOL®

Dermacol is a unique stoma collar that provides a leakproof barrier around the base of your stoma, preventing output from coming into contact with the skin. Unlike any other product, the Dermacol stoma collar forms a physical barrier to prevent output from coming into contact with the delicate skin around your stoma, protecting you from discomfort and giving you extra peace of mind. It's suitable for use with all types of appliances.

Skin-friendly hydrocolloid



Ideal size and height

Flexible collar




DERMACOL®		
STOMA SIZE	PRODUCT CODE	PACK SIZE
20mm – 23mm	DC23	30
24mm – 26mm	DC26	30
27mm – 29mm	DC29	30
30mm – 32mm	DC32	30
33mm – 35mm	DC35	30
36mm – 38mm	DC38	30

SECUPLAST® MOULDABLE SEALS®

New SecuPlast Mouldable Seals help prevent sore skin by reducing the leaks that are the main cause of it.

Simply mould SecuPlast Mouldable Seals to the exact size and shape of your stoma and the exact contours of your abdomen to create a customised, comfortable and leakproof seal.

- No need to warm
- Suitable for any make or style of pouch and all stoma types
- Very easy to mould and shape to fit; then place in position
- Very easy to rejoin if it breaks
- Easy to break into smaller pieces, if required
- Skin-friendly, with no sting, even on sore skin
- Customised comfort – SecuPlast Mouldable Seals keep the shape you create for secure and lasting protection from leaks




AVAILABLE IN

Large

Standard

Thin



SECUPLAST® MOULDABLE SEALS

SIZE	DIAMETER	THICKNESS	ORDER CODE	PACK SIZE
THIN	50mm	3mm	SMST	30
STANDARD	50mm	4.2mm	SMSS	30
LARGE	100mm	3mm	SMSL	10

