

# STOMA SOLUTIONS

A range of easy-to-use products ensuring comfort and peace of mind



All the products featured in this brochure are available through Salts Medilink® or your usual supplier of ostomy products.

For more information and advice about Stoma Solutions or the Salts Medilink® 24\* Hour Home Delivery Service, please contact:

**FREEPHONE 0800 626388** (UK) or **FREEPHONE 1-800 408508** (Ireland)  
or log on to our website at [www.salts.co.uk/medilink](http://www.salts.co.uk/medilink)



100% designed and manufactured in the UK

**For free samples please call**

**FREEPHONE 0800 626388** (UK) or

**FREEPHONE 1-800 408508** (Ireland)

or log on to our website at

[www.saltsstomacare.co.uk](http://www.saltsstomacare.co.uk)

**SALTS**  
HEALTHCARE

© Registered trademarks of Salts Healthcare Ltd. © Salts Healthcare Ltd 2010. Products and trade marks of Salts Healthcare Limited are protected by UK and foreign patents, registered designs and trade marks. For further details, please visit [www.salts.co.uk](http://www.salts.co.uk) \*If required, except Highlands and Islands and subject to availability. RM763888 04/10/3K



Validating Salts research  
into healthy stoma skin

**SALTS**  
HEALTHCARE

Excellence in stoma care

# STOMA SOLUTIONS

## The solution for any problem

Our easy-to-use range has been specifically developed to alleviate the problems faced by ostomists, to give you security, confidence and peace of mind.

The problems that our accessories are designed to address are: **LEAKS:** the most common problem that invariably leads to others. **SORE SKIN:** due to leakage and constant changing of the appliance. **ADHESION:** lack of adhesion caused by uneven skin around the stoma. **ODOUR:** caused by bacteria in the output of the stoma. **RESIDUE:** surplus adhesive left on the skin. **SKIN PROTECTION:** for maximum comfort and healing. We also provide a number of **OTHER PRODUCTS.** Your Stoma Nurse will be able to recommend the appropriate product and should always be consulted regarding any problems that you may have.

### FREE SAMPLES

Try Stoma Solutions and see the difference. **For free samples please call FREEPHONE 0800 626388 (UK) or FREEPHONE 1-800 408508 (Ireland) or log on to our website at [www.saltsstomacare.co.uk](http://www.saltsstomacare.co.uk)**

**Stoma Solutions** can be ordered easily from Salts Medilink® by post, phone or online. Alternatively, order from your usual supplier of ostomy products.

For more information and advice about Stoma Solutions or the Salts Medilink® 24\* Hour Home Delivery Service, please contact: **FREEPHONE 0800 626388 (UK) or FREEPHONE 1-800 408508 (Ireland) or log on to our website at [www.salts.co.uk](http://www.salts.co.uk)**



# Simplicity® STOMA SUPPORT WEAR

LEAKS

SORE SKIN

ADHESION


ODOUR

RESIDUE

SKIN PROTECTION

OTHER PRODUCTS

Available in both black and white, our ladies' briefs and unisex boxer are made from a **two-way stretch fabric** that allows the body to move freely, while ribbing on the front provides a flush finish.



- Ribbing, for secure pouch support and flush finish
- A higher waistband to cover pouch
- Deep waistband, for extra comfort
- Virtually invisible under most clothing
- Soft and comfortable against the skin
- Lightweight and breathable antibacterial fabric



### SIMPLICITY® UNISEX BOXER

SIZE	COLOUR	ORDER CODES
X Small (22-27")	Black	BOXBXS
X Small (22-27")	White	BOXWXS
Small/Med (28/30/32")	Black	BOXBSM
Small/Med (28/30/32")	White	BOXWSM
Med/Large (34/36/38")	Black	BOXBML
Med/Large (34/36/38")	White	BOXWML
Large/X Large (40/42/44")	Black	BOXBLXL
Large/X Large (40/42/44")	White	BOXWLXL



### SIMPLICITY® LADIES BRIEF

SIZE	COLOUR	ORDER CODES
X Small (4/6/8)	Black	BRFBXS
X Small (4/6/8)	White	BRFWXS
Small/Med (8/10/12)	Black	BRFBSM
Small/Med (8/10/12)	White	BRFWSM
Med/Large (14/16/18)	Black	BRFBML
Med/Large (14/16/18)	White	BRFWML
Large/X Large (20/22/24)	Black	BRFBLXL
Large/X Large (20/22/24)	White	BRFWLXL



# PERI-PREP® SENSITIVE


**Dermacol®**  
stoma collar FOR COMPLETE SECURITY

LEAKS SORE SKIN ADHESION ODOUR RESIDUE **SKIN PROTECTION** OTHER PRODUCTS

Leakage and the trauma caused by the constant changing of a pouch can lead to sore skin. Peri-Prep Sensitive is specifically designed to leave a gentle protective film on the skin immediately around the stoma.



**New technology:** Peri-Prep Sensitive's formulation, developed in conjunction with Aston University, provides a durable, leak-resistant barrier when dried, but is pH neutral to be gentle on the skin. The composition contains 'star shaped' polymers that allow greater flexibility of the film and improved pouch adhesion.



- Prevents sore skin following frequent pouch changes
- Prevents perspiration penetrating pouch adhesive
- Does not need to be removed on every pouch change
- Protects skin against leakages


If required, Peri-Prep Sensitive can be built up into layers by repeated application to provide a long-lasting barrier, for extra security. Peri-Prep Sensitive is gentle on the skin and provides maximum protection – without stinging.

PERI-PREP® SENSITIVE		
SIZE	QUANTITY	ORDER CODE
–	30 sachets	PPS1

LEAKS SORE SKIN **ADHESION** ODOUR RESIDUE **SKIN PROTECTION** OTHER PRODUCTS

Dermacol is a unique stoma collar that provides a leakproof barrier around the base of your stoma, preventing output from coming into contact with the skin.

Unlike any other product, the Dermacol stoma collar forms a physical barrier to prevent output from coming into contact with the delicate skin around your stoma, protecting you from discomfort and giving you extra peace of mind. It's suitable for use with all types of appliances.



DERMACOL®		
STOMA SIZE	PRODUCT CODE	PACK SIZE
20mm – 23mm	DC23	30
24mm – 26mm	DC26	30
27mm – 29mm	DC29	30
30mm – 32mm	DC32	30
33mm – 35mm	DC35	30
36mm – 38mm	DC38	30

LEAKS

SORE SKIN

ADHESION

ODOUR

RESIDUE

SKIN  
PROTECTIONOTHER  
PRODUCTS

LEAKS

SORE SKIN

ADHESION

ODOUR

RESIDUE

SKIN  
PROTECTIONOTHER  
PRODUCTS

You may occasionally experience adhesion difficulties with your pouch. For extra security, the simple and effective solution is SecuPlast Hydro. Forget sticky and uncomfortable tape – SecuPlast Hydro strips are thin, skin friendly and very easy to apply.

- Extra security with any type or make of pouch
- Easy to use – allows precise positioning and less wastage, so you can use as many strips as you require
- Skin-friendly – unlike microporous tapes
- Non-allergenic and can be used by patients who cannot use microporous tapes
- One size fits all


All the security you need, wherever you need it. Don't just take our word for it. After extensive product testing:

- 96% of our patients said SecuPlast Hydro was comfortable
- 74% said the size was about right
- 62% found SecuPlast Hydro better than their usual product



New SecuPlast Mouldable Seals help prevent sore skin by reducing the leaks that are the main cause of it.

Simply mould SecuPlast Mouldable Seals to the exact size and shape of your stoma and the exact contours of your abdomen to create a customised, comfortable and leakproof seal.



- No need to warm
- Suitable for any make or style of pouch and all stoma types
- Very easy to mould and shape to fit; then place in position
- Very easy to rejoin if it breaks
- Easy to break into smaller pieces, if required
- Skin-friendly, with no sting, even on sore skin
- Customised comfort – SecuPlast Mouldable Seals keep the shape you create for secure and lasting protection from leaks



## SECUPLAST<sup>®</sup> HYDRO

SIZE	QUANTITY	ORDER CODE
One size fits all	30	SPH1

## SECUPLAST<sup>®</sup> MOULDABLE SEALS

SIZE	DIAMETER	THICKNESS	ORDER CODE	PACK SIZE
THIN	50mm	3mm	SMST	30
STANDARD	50mm	4.2mm	SMSS	30
LARGE	100mm	3mm	SMSL	10

# WIPEAWAY ADHESIVE REMOVER

# FRESHAIRE ODOUR ELIMINATING SPRAY

LEAKS SORE SKIN ADHESION ODOUR RESIDUE SKIN PROTECTION OTHER PRODUCTS LEAKS SORE SKIN ADHESION ODOUR RESIDUE SKIN PROTECTION OTHER PRODUCTS

Tried and trusted WipeAway adhesive remover. For years, our popular WipeAway sachets have been used to remove pouches and adhesive residue and to help keep the skin around the stoma healthy.



The natural ingredients help to soothe skin trauma associated with adhesive removal, leaving a hint of fresh lavender.

- Stress-free removal of pouches and adhesive residue
- Reduces trauma associated with adhesive removal
- Skin-friendly: does not contain alcohol

**Do not wipe directly onto stoma.**  
**If necessary, pat dry.**

Unlike many other air fresheners and deodorisers, FreshAire doesn't just mask unpleasant smells with perfume, but attacks and destroys the odour molecules completely, so the odour can't come back.

- Eliminates odour – it really does work!
- Odours are eliminated at a molecular level
- Not just another masking or neutralising spray
- Non-toxic, non-hazardous and non-flammable
- Environmentally friendly
- Non-allergenic

FreshAire is effective and easy to use. During pouch changes, FreshAire is simply sprayed once or twice into the air around the pouch.

- Lightweight and a handy size
- Can be used safely in confined spaces
- Will not harm skin surfaces or fabrics
- A hint of peppermint



WIPEAWAY ADHESIVE REMOVER			
DESCRIPTION	SIZE	QUANTITY	ORDER CODE
Wipes	–	30 sachets	WA1

FRESHAIRE ODOUR ELIMINATING SPRAY		
SIZE	PACK SIZE	ORDER CODE
50ml	1	FA1

# STOMA PASTE

# PROTECTIVE WAFERS

LEAKS

SORE SKIN

ADHESION

ODOUR

RESIDUE

SKIN PROTECTION

OTHER PRODUCTS

LEAKS

SORE SKIN

ADHESION

ODOUR

RESIDUE

SKIN PROTECTION

OTHER PRODUCTS

Leakage can often occur when the skin immediately around the stoma is uneven. These gaps prevent a perfect seal and cause leaks that can eventually lead to sore skin.

- Creates a perfect seal, for maximum adhesion and increased security
- Improves pouch fit around the stoma
- 30 seconds drying time
- Clean and easy to use

Stoma Paste has been specifically designed to be used as an effective filler where there is scar tissue or gaps around the stoma. Stoma Paste is easy to use and, once applied, can be spread into position, providing an ideal surface for appliances to adhere to.

With a drying time of just 30 seconds, Stoma Paste is not messy to use and will not stick to fingers, unlike many of the other pastes available.



An allergic or undesirable reaction to an adhesive used on a particular appliance may be the cause of sore skin. Indeed, the adhesive itself may not even be strong enough to fully secure the appliance.



- Non-irritating and suitable for sensitive skin
- Thinner than other leading brands, for maximum discretion
- Easy to cut and simple to apply
- Promotes healing when used on broken skin
- Flexible and soft, for added comfort

Protective Wafers are specifically designed to help prevent soreness and increase adhesion. Salts' hydrocolloid wafers provide an ideal environment to promote healing of sore or broken skin – and then help prevent the condition from reoccurring.

STOMA PASTE		
SIZE	QUANTITY	ORDER CODE
60g	1	SP60

PROTECTIVE WAFERS		
SIZE	QUANTITY	ORDER CODE
10cm x 10cm	10	PW1010
15cm x 15cm	5	PW1515

# NO-ROMA

# SALTAIR LOTION AND SOAP

LEAKS

SORE SKIN

ADHESION

ODOUR

RESIDUE

SKIN PROTECTION

OTHER PRODUCTS

LEAKS

SORE SKIN

ADHESION

ODOUR

RESIDUE

SKIN PROTECTION

OTHER PRODUCTS

Unpleasant odour escaping from an appliance can be very embarrassing. No-Roma, the easy-to-use pouch deodorant, is specifically designed to prevent this situation from occurring.

- Kills the odour-producing bacteria
- Available in two sizes – pocket size is ideal to carry around
- Completely odourless
- Lasting results

Unlike other products, No-Roma works by actually killing the bacteria that produce the odour, rather than just simply masking it with perfume. It is also very simple to use. When changing an appliance, a drop of No-Roma is placed inside the pouch (not on the stoma itself) to prevent the build-up of odour.



No-Roma can be used with any appliance, and is available in a small bottle, making it ideal to discreetly carry around.

## SALTAIR LOTION ●

Skin irritation or soreness can be uncomfortable and distressing. SaltAir Lotion helps to soothe sore or irritated skin – and because it contains no alcohol, it can be used on sore skin without stinging.

- Easy to use
- Helps soothe sore or irritated skin
- Alcohol-free, to eliminate stinging
- Non-greasy
- Free of charge with orders placed through Salts Medilink®

SaltAir Lotion is quick and easy to use, by simply applying to a dry wipe and dabbing onto the affected area. It is not greasy and helps the pouch to adhere better to the skin when fully dry. SaltAir Lotion can be used on skin around any type of stoma – ileostomy, colostomy or urostomy – under any type of appliance. *SaltAir Lotion is only available through Salts Medilink®*



## SALTAIR SOAP ●

On some occasions, warm water alone may not be enough to remove adhesive residue left on the skin. In other cases, the skin may also be slightly inflamed, sore or even broken.

- Quickly removes adhesive residue
- Cleans skin and washes off easily
- Non-perfumed
- Easy to use

SaltAir Soap is a soft, liquid soap designed specifically for ostomists for situations where gentle cleansing is required. SaltAir Soap is non-perfumed to reduce the risk of irritating the skin. It is also quick and easy to use, with just a few drops added to water before wiping the skin.

SaltAir Soap quickly removes excess adhesive without the need for rubbing and rinses away without leaving a residue.



### NO-ROMA

SIZE	QUANTITY	ORDER CODE
28ml	1	833021
227ml	1	833056

### SALTAIR LOTION AND SOAP

DESCRIPTION	SIZE	QUANTITY	ORDER CODE
Lotion	28ml	1	833006SAMP
Soap	110ml	1	833007

