

UROSTOMY PRODUCTS



All the products featured in this brochure are available through Salts Medilink® or your usual supplier of ostomy products.

For more information and advice about products featured in this brochure, Stoma Solutions or the Salts Medilink® 24* Hour Home Delivery Service, please contact:

FREEPHONE 0800 626388 (UK) or **FREEPHONE 1-800 408508** (Ireland)
or log on to our website at www.salts.co.uk/medilink



100% designed and manufactured in the UK

For free samples please call

FREEPHONE 0800 626388 (UK) or

FREEPHONE 1-800 408508 (Ireland)

or log on to our website at

www.saltsstomacare.co.uk

SALTS
HEALTHCARE

© Registered trademarks of Salts Healthcare Ltd. © Salts Healthcare Ltd 2010. Products and trade marks of Salts Healthcare Limited are protected by UK and foreign patents, registered designs and trade marks. For further details, please visit www.salts.co.uk *If required, except Highlands and Islands and subject to availability. RM763887 04/10/3K



A complete guide to our comprehensive range of urostomy products for the one- and two-piece user



Validating Salts research
into healthy stoma skin

SALTS
HEALTHCARE

Excellence in stoma care

Our top-selling range of urostomy pouches for the one-piece user, incorporating a soft and comfortable beige material on both sides.



EASY TO USE

A simple, single-release clear film with a peel-off tab for quick and easy application, and a split comfort material, for easy positioning.

SKIN-FRIENDLY

A proven, non-allergenic, soft and flexible hydrocolloid wafer, for even the most sensitive skin.

EXTRA COMFORT

A soft and comfortable beige material on both sides.

SECURITY

A generous, diamond-shaped hydrocolloid wafer, for added security and peace of mind.

DISCREET

Low profile, for better cosmetic appearance.



TRANSPARENT/OVERLAP

STANDARD POUCHES

SIZE	STARTER HOLE	PRE-CUT	ORDER CODES			PACK SIZE
			Transparent	Trans/Overlap	Opaque	
13mm	•	–	–	UNW13	–	10
13mm*	•	–	–	UTO13	–	10
25mm	–	•	–	UNW25	–	10
32mm	–	•	–	UNW32	–	10
38mm	–	•	–	UNW38	–	10

SMALL POUCH

SIZE	STARTER HOLE	PRE-CUT	ORDER CODES			PACK SIZE
			Transparent	Trans/Overlap	Opaque	
13mm	•	–	–	USNW13	–	10

* Oval wafer pouches can be cut up to 110mm wide and 70mm deep.



* Oval wafer.



Positive action, wide-bore on/off tap with integral drip cap.



Small pouch.

A convex pouch with even more flexibility and softness. Confidence Convex Supersoft incorporates a soft, flexible wafer with a gentle profile – for increased comfort.

1
PIECE

SOFT

More flexible, softer convex flange for added comfort.

COMFORTABLE

Shallower profile fits more comfortably against the skin.

FLEXIBLE

New, flexible foam backing conforms better to body shape.

SECURE

New hydrocolloid wafer with 'high tack' adhesive affords extra security.

CUSHIONED

Soft padding acts as a cushion between skin and the on/off tap.

EASY-TO-USE

Positive action, wide-bore on/off tap with integral drip cap.



STANDARD POUCHES				
SIZE	STARTER HOLE	PRE-CUT	ORDER CODES	PACK SIZE
13–25mm	•	–	CUSS1325	10
13–38mm	•	–	CUSS1338	10
13–52mm	•	–	CUSS1352	10
21mm	–	•	CUSS21	10
25mm	–	•	CUSS25	10
28mm	–	•	CUSS28	10
32mm	–	•	CUSS32	10
35mm	–	•	CUSS35	10
38mm	–	•	CUSS38	10



Overlapping fabric for easy positioning.



Positive action, wide-bore on/off tap with integral drip cap.



This product, with integral convexity, should only be used after assessment by a medical professional.

A range of urostomy pouches for the two-piece user, incorporating a soft, silky cover of breathable material.

UNIQUE CONNECTION

Unique flange connection, with audible 'click' to ensure that the pouch is fitted correctly for extra security.

ROTATING POUCH

Pouches can be rotated for maximum comfort.



BODYSIDE WAFERS

For use in conjunction with Second Nature two-piece pouches.



Overlapping fabric for easy positioning.



Positive action, wide-bore on/off tap with integral drip cap.

STANDARD POUCHES				BODYSIDE WAFERS		
SIZE	TO FIT STOMA SIZE UP TO	ORDER CODES Trans/Overlap	PACK SIZE	WAFER SIZE	ORDER CODES	PACK SIZE
32mm	26mm	2U32	10	110 x 100mm	2FL32	5
38mm	32mm	2U38	10	110 x 100mm	2FL38	5
45mm	39mm	2U45	10	110 x 100mm	2FL45	5
57mm	51mm	2U57	10	134 x 124mm	2FL57	5

A range of urostomy pouches for the one-piece user, incorporating a soft, silky cover of breathable material.

EASY TO USE

A simple, single-release clear film with a peel-off tab, for quick and easy application and a split comfort material, for easy positioning.

SECURITY

A generous, diamond-shaped wafer, for added security and peace of mind.

DISCREET

Low profile, for a better cosmetic appearance.



SKIN-FRIENDLY

A proven, non-allergenic, soft and flexible hydrocolloid wafer, for even the most sensitive skin.

COMFORT

A soft, silky cover of breathable material, for extra comfort. The water-resistant material can be patted dry after bathing or showering.

UNIQUE DESIGN FEATURE

Positive action, wide-bore open/close tap, for easy draining, with integral drip cap.

STANDARD POUCHES						
SIZE	STARTER HOLE	PRE-CUT	ORDER CODES		PACK SIZE	
			Transparent	Trans/Overlap Opaque		
13mm	•	–	–	U13	10	
25mm	–	•	–	U25	10	
32mm	–	•	–	U32	10	
38mm	–	•	–	U38	10	

STOMA SOLUTIONS

The solution for any problem. Our easy-to-use range has been specifically developed to alleviate the problems faced by ostomists, to provide security, confidence and peace of mind. Please see the Stoma Solutions brochure for the full range of accessories.



NIGHT DRAINAGE BAG

- Longer 1,200mm tube allows greater movement
- Flexible material made from anti-kink material
- Quick and easy snap-on connection
- Unique swivel action connection allows 360° rotation
- Large 2,000ml capacity



SECUPLAST® HYDRO

Skin-friendly hydrocolloid strips to provide extra security around the edge of the flange.



PERI-PREP® SENSITIVE

Leakage and the trauma caused by the constant changing of a pouch can lead to sore skin. Peri-Prep Sensitive is specifically designed to leave a gentle protective film on the skin immediately around the stoma.




WIPEAWAY

WipeAway sachets are used to remove pouches and adhesive residue and to help keep the skin around the stoma healthy. The natural ingredients help to soothe skin trauma associated with adhesive removal, leaving a hint of fresh lavender.



SIMPLICITY STOMA SUPPORT WEAR

Available in both black and white, our ladies' briefs and unisex boxer are made from a two-way stretch fabric that allows the body to move freely, while ribbing on the front provides a flush finish.



NIGHT DRAINAGE BAG & ADAPTOR			SECUPLAST® HYDRO		
SIZE	ORDER CODE	PACK SIZE	SIZE	ORDER CODE	PACK SIZE
2,000ml	ZL0400	10	One size fits all	SPH1	30
Adaptor	NDA6	6			
PERI-PREP® SENSITIVE			WIPEAWAY		
SIZE	ORDER CODE	PACK SIZE	SIZE	ORDER CODE	PACK SIZE
–	PPS1	30 sachets	–	WA1	30 sachets

SIMPLICITY® UNISEX BOXER			SIMPLICITY® LADIES BRIEF		
SIZE	COLOUR	ORDER CODES	SIZE	COLOUR	ORDER CODES
X Small (22–27")	Black	BOXBXS	X Small (4/6/8)	Black	BRFBXS
X Small (22–27")	White	BOXWXS	X Small (4/6/8)	White	BRFWXS
Small/Med (28/30/32")	Black	BOXBSM	Small/Med (8/10/12)	Black	BRFBSM
Small/Med (28/30/32")	White	BOXWSM	Small/Med (8/10/12)	White	BRFWSM
Med/Large (34/36/38")	Black	BOXBML	Med/Large (14/16/18)	Black	BRFBML
Med/Large (34/36/38")	White	BOXWML	Med/Large (14/16/18)	White	BRFWML
Large/X Large (40/42/44")	Black	BOXBLXL	Large/X Large (20/22/24)	Black	BRFBLXL
Large/X Large (40/42/44")	White	BOXWLXL	Large/X Large (20/22/24)	White	BRFWLXL

DERMACOL®

Dermacol is a unique stoma collar that provides a leakproof barrier around the base of your stoma, preventing output from coming into contact with the skin. Unlike any other product, the Dermacol stoma collar forms a physical barrier to prevent output from coming into contact with the delicate skin around your stoma, protecting you from discomfort and giving you extra peace of mind. It's suitable for use with all types of appliances.



Ideal size and height

Flexible collar

Skin-friendly hydrocolloid




DERMACOL®		
STOMA SIZE	PRODUCT CODE	PACK SIZE
20mm – 23mm	DC23	30
24mm – 26mm	DC26	30
27mm – 29mm	DC29	30
30mm – 32mm	DC32	30
33mm – 35mm	DC35	30
36mm – 38mm	DC38	30

SECUPLAST® MOULDABLE SEALS®

New SecuPlast Mouldable Seals help prevent sore skin by reducing the leaks that are the main cause of it.

Simply mould SecuPlast Mouldable Seals to the exact size and shape of your stoma and the exact contours of your abdomen to create a customised, comfortable and leakproof seal.



- No need to warm
- Suitable for any make or style of pouch and all stoma types
- Very easy to mould and shape to fit; then place in position
- Very easy to rejoin if it breaks
- Easy to break into smaller pieces, if required
- Skin-friendly, with no sting, even on sore skin
- Customised comfort – SecuPlast Mouldable Seals keep the shape you create for secure and lasting protection from leaks



SECUPLAST® MOULDABLE SEALS				
SIZE	DIAMETER	THICKNESS	ORDER CODE	PACK SIZE
THIN	50mm	3mm	SMST	30
STANDARD	50mm	4.2mm	SMSS	30
LARGE	100mm	3mm	SMSL	10